

3M Science.
Applied to Life.™





3M™ RelyX™ Universal Resin Cement &
3M™ Scotchbond™ Universal Plus Adhesive

The two that can do it all.

3M™ RelyX™ Universal Resin Cement works both as a standalone, self-adhesive cement – and as an adhesive cement together with 3M™ Scotchbond™ Universal Plus Adhesive. The adhesive serves as a universal primer for all materials and covers all direct and indirect bonding indications.

Eliminate the hassle and confusion of multiple resin cements, primers and adhesives and simplify your direct and indirect restorative workflows. Fewer products to keep in stock and replenish. Saves space, time and money.




-  **First radiopaque all-in-one universal adhesive**
-  **Further boosts bond strength of RelyX Universal Resin Cement**
-  **Fully aligned system: adhesive cured by cement, no light-cure needed**
-  **Universal primer for all restorative materials**



Game-changing syringe. Proprietary 3M design.

-  **Micro Mixing Tip providing 80% less cement waste***
-  **15 applications from 3.4g****
-  **More hygienic, easier to clean***** 
self-sealing
-  **More ergonomic shape with 50% less plastic waste***
-  **Thin, long elongation tip for easy root canal application** 
-  **Easy excess clean-up** 
Image courtesy of Dr. Jan Gueth, Germany
-  **Superior self-adhesive bond strength to dentin in self- and light-cure**


Shear bond strength (MPa) to dentin after artificial aging



R. Afutu, M. Abreu, G. Kugel, Tufts University, J. Dent. Res. Vol 98A, No 3629, 2019

* per application compared to currently available standard automix systems
** on average *** compared to currently used automix syringes

High esthetics with four color-stable, fluorescent shades



Ordering information

3M™ RelyX™ Universal Resin Cement

Item no. Product information

Intro Kit

56968 4x 3M™ RelyX™ Universal Resin Cement syringes shades TR, A1, A3O, WO (3.4g each), 60x 3M™ RelyX™ Universal Micro Mixing Tips, 20 Elongation Tips, 1 vial 3M™ Scotchbond™ Universal Plus Adhesive (5 ml), 1x 3M™ Scotchbond™ Universal Etchant (3 ml), 25 blue etchant tips, 50 orange microbrushes, step-by-step guides

Trial Kits

Each trial kit contains: 1x 3M™ RelyX™ Universal Resin Cement syringe (3.4g), 15x 3M™ RelyX™ Universal Micro Mixing Tips, 5 Elongation Tips, 1 vial 3M™ Scotchbond™ Universal Plus Adhesive (1.5 ml), 1x 3M™ Scotchbond™ Universal Etchant (3 ml), 25 blue etchant tips, 50 orange microbrushes, step-by-step guides

56969 Shade TR 56970 Shade A1

Refills

Each refill contains: 1x 3M™ RelyX™ Universal Resin Cement syringe (3.4g), 15x 3M™ RelyX™ Universal Micro Mixing Tips, 5 Elongation Tips

56971 Shade TR 56973 Shade A3O
56972 Shade A1 56974 Shade WO

Value Packs

Each value pack contains: 3x 3M™ RelyX™ Universal Resin Cement syringes (3.4g each), 45x 3M™ RelyX™ Universal Micro Mixing Tips, 15 Elongation Tips

56977 Shade TR 56978 Shade A1

Mixing Tips

56975 30x 3M™ RelyX™ Universal Micro Mixing Tips
56976 15x 3M™ RelyX™ Universal Micro Mixing Tips with Elongation Tips

3M™ Scotchbond™ Universal Plus Adhesive

Vial

41293 Intro Kit: 1x 3M™ Scotchbond™ Universal Plus Adhesive vial (5 ml), 50 disposable brushes, 1x 3M™ Scotchbond™ Universal Etchant (3 ml), 25 dispensing tips

Refill: 3M™ Scotchbond™ Universal Plus Adhesive vial (5 ml)

41294 1x 5 ml vial 41295 3x 5 ml vial

Unit Dose

3M™ Scotchbond™ Universal Plus Adhesive Unit Dose, 1x 3M™ Scotchbond™ Universal Etchant (3 ml), 25 dispensing tips

41297 Intro Kit: 50x L-Pop 41299 200x L-Pop
41298 100x L-Pop 41304 400x L-Pop

3M™ Scotchbond™ Universal Etchant Etching Gel

41263 2x 3M™ Scotchbond™ Universal Etchant (3 ml each), 50 dispensing tips



3M™ RelyX™ Universal Resin Cement Intro Kit (56968)



3M™ RelyX™ Universal Resin Cement Trial Kit TR (56969)



3M™ RelyX™ Universal Resin Cement Refill TR (56971)



3M™ Scotchbond™ Universal Plus Adhesive Refill (41294)



3M™ Scotchbond™ Universal Plus Adhesive Unit Dose:

For efficient hygiene management. Reduces the risk of cross-contamination.

3M™ Scotchbond™ Universal Plus Adhesive Unit Dose 100 L-Pop (41298)

3M.com/RelyXUniversal

3M.com/ScotchbondPlus

