

CEMENTS & LINERS



Cement selection made easy.

	3M™ RelyX™ Universal Resin Cement (see statement below)	3M™ RelyX™ Ultimate Adhesive Resin Cement	3M™ RelyX™ Unicem/Unicem 2 Self-Adhesive Resin Cement	3M™ RelyX™ Luting Plus Resin Modified Glass Ionomer Cement	3M™ RelyX™ Veneer Cement
Metal Based	(on pg. 88)	(on pg. 86)	(on pg. 86)	(on pg. 79)	(on pg. 88)
Inlays/Onlays	++	+	++	+	-
Crowns/Bridges	++	+	++	++	-
Endodontic Posts	++	+	++	++	-
Maryland Bridges	++	++	+	-	-
On Implant Abutments	++	+	++	++	-
Glass Ceramics including Lithium Disilicate (e.g. IPS e.max®, VITABLOCS® Mark II, Celtra® Duo, IPS Empress®)					
Inlays/Onlays/Table Tops	++	++	+	-	-
Crowns/Bridges	++	++	++	-	-
Veneers	+	+	-	-	++
On Implant Abutments	++	+	+	-	-
Zirconia and other Oxide Ceramics (e.g. 3M™ Lava™ Zirconia, BruxZir® Zirconia)					
Inlays/Onlays	++	++	++	+	-
Crowns/Bridges	++	+	++	+	-
Endodontic Posts	++	+	++	+	-
Maryland Bridges	++	++	+	-	-
On Implant Abutments	++	+	++	++	-
Resin Based Composites/ Hybrids (e.g. Cerasmart, Enamic, 3M™ Lava™ Ultimate CAD/CAM Restorative); Fiber Reinforced Composite Posts (e.g. RelyX™ Fiber Posts)					
Inlays/Onlays/Table Tops	++	++	+	-	-
Crowns	++	++	++	-	-
Endodontic Posts	++	+	++	-	-
Veneers	+	+	-	-	++
On Implant Abutments	++	++	++	-	-
++ Highly recommended* + Recommended - Not indicated * Either better performance for this indication or easier handling with equal performance.					

Note: The cement selection recommendations in this guide list the key indications for a particular cement, and are intended as a general guide only. 3M RelyX cements are highly versatile products that can be used for multiple indications. Your selection is dependent on individual patient circumstances and your preferred procedural techniques.

Depending on the indication, 3M™ RelyX™ Universal Resin Cement can be used as a self-adhesive resin cement or as an adhesive resin cement when combined with 3M™ Scotchbond™ Universal Plus Adhesive. Please refer to 3M™ RelyX™ Universal Resin Cement instructions for use for details.

Learn more at 3M.com/dentalcements

