

Alpen Conversion List

Important Information	EN
-----------------------	----

We are delighted to inform our valued clients that COLTENE has made significant improvements regarding the packaging of our Alpen Diamonds:

- User-friendly blister packaging format
- Securely sealed, preserving integrity of the product
- Equipped with an easy peel feature

Additionally, we have made certain changes in ring colors and grit sizes to align with the ISO 7711 standard (see table below):

- The black "super coarse" (MLX) has been reclassified as green "coarse" (C)
- The green "coarse" (ML) has been reclassified as blue "medium" (M)

Furthermore, there have been modifications to the letters in the names of our burs:

- The black MLX instruments have been transformed to green C
- The blue ML to blue M
- The yellow XF has been replaced by yellow EF



There have been changes to the item numbers. For your reference, we have included a conversion table that lists the old and new item numbers. This conversion list only includes new items that have an equivalent old item. However, the new portfolio also includes completely new items not shown in this list. You can access the PDF version of this table through the provided link: <https://online.flippingbook.com/view/855353489/> or by scanning the QR code.

We appreciate your understanding as we strive to ensure precision and quality and provide the best possible solutions for you.

Your COLTENE Team

Old color ring - index (grit type)	Old mean grit value	New color ring - index (grit type)	New grit value range and mean
Black - SC (super coarse)	140 µm	Green - C (coarse)	107-181 µm; Mean 151 µm
Green - C (coarse)	120 µm	Blue - M (medium)	64-126 µm; Mean 107 µm
Grey - M (medium)	105 µm	Blue - M (medium)	64-126 µm; Mean 107 µm
Red - F (fine)	40 µm	Red - F (fine)	27-76 µm; Mean 46 µm
Yellow - EF (extra fine)	30 µm	Yellow - EF (extra fine)	10-36 µm; Mean 25 µm
White - UF (ultra fine)	20 µm	not available anymore	4-14 µm; Mean 8 µm
Black - TSX (Turbo Cut)	140µm	Green - C (coarse) or Blue - M (medium)	does not exist according to ISO

REF new	Description new	REF old	Description old
60031753	G368-314-016-03.5-EF, Alpen Diamonds	R368EF016FG	G368-314-016-3.5-XF
60031755	G368-314-016-03.5-F, Alpen Diamonds	R368F016FG	G368-314-016-3.5-F
60031757	G368-314-016-03.5-M, Alpen Diamonds	R368C016FG	G368-314-016-3.5-ML
60031759	G368-314-020-05.0-C, Alpen Diamonds	R368SC020FG	G368-314-020-5-MLX
60031761	G368-314-020-05.0-F, Alpen Diamonds	R368EF020FG	G368-314-020-5-XF
60031763	G368-314-020-05.0-F, Alpen Diamonds	R368F020FG	G368-314-020-5-F
60031767	G368-314-020-05.0-M, Alpen Diamonds	R368C020FG	G368-314-020-5-ML
60031769	G368-314-023-05.0-C, Alpen Diamonds	R368SC023FG	G368-314-023-5-MLX
60031771	G368-314-023-05.0-EF, Alpen Diamonds	R368EF023FG	G368-314-023-5-XF
60031773	G368-314-023-05.0-F, Alpen Diamonds	R368F023FG	G368-314-023-5-F
60031777	G368-314-023-05.0-M, Alpen Diamonds	R368C023FG	G368-314-023-5-ML
60031781	G379-314-016-03.5-EF, Alpen Diamonds	R379FE016FG	G379-314-016-3.5-XF
60031783	G379-314-016-03.5-F, Alpen Diamonds	R379F016FG	G379-314-016-3.5-F
60031785	G379-314-016-03.5-M, Alpen Diamonds	R379C016FG	G379-314-016-3.5-ML
60031787	G379-314-018-03.5-EF, Alpen Diamonds	R379FE018FG	G379-314-018-3.5-XF
60031789	G379-314-018-03.5-F, Alpen Diamonds	R379F018FG	G379-314-018-3.5-F
60031793	G379-314-018-03.5-M, Alpen Diamonds	R379C018FG	G379-314-018-3.5-ML
60031795	G379-314-023-05.0-C, Alpen Diamonds	R379SC023FG	G379-314-023-5-MLX
60031797	G379-314-023-05.0-EF, Alpen Diamonds	R379FE023FG	G379-314-023-5-XF
60031801	G379-314-023-05.0-F, Alpen Diamonds	R379F023FG	G379-314-023-5-F
60031807	G379-314-023-05.0-M, Alpen Diamonds	R379C023FG	G379-314-023-5-ML
60031811	G379T-314-023-05.0-C, Alpen Diamonds	R379TC023FG	G379-314-023-5-TSX
60031813	G379Z-314-023-05.0-M, Alpen Diamonds	ZREX379M023FG	ZR379-314-023-5-M
60031815	G390-314-014-03.0-EF, Alpen Diamonds	R390FE014FG	G390-314-014-3-XF
60031817	G390-314-014-03.0-F, Alpen Diamonds	R390F014FG	G390-314-014-3-F
60031819	G392-314-016-06.0-EF, Alpen Diamonds	R392FE016FG	G392-314-016-6-XF
60031821	G392-314-016-06.0-F, Alpen Diamonds	R392F016FG	G392-314-016-6-F

REF new	Description new	REF old	Description old
60031825	G392-314-016-06.0-M, Alpen Diamonds	R392C016FG	G392-314-016-6-ML
60031827	G392-314-016-08.0-F, Alpen Diamonds	R3928F016FG	G392-314-016-8-F
60031829	G392-314-016-08.0-M, Alpen Diamonds	R3928C016FG	G392-314-016-8-ML
60031831	G392-314-018-08.0-EF, Alpen Diamonds	R392EF018FG	G392-314-018-8-XF
60031833	G392-314-018-08.0-F, Alpen Diamonds	R392F018FG	G392-314-018-8-F
60031835	G392-314-018-08.0-M, Alpen Diamonds	R392C018FG	G392-314-018-8-ML
60031837	G392-314-021-10.0-F, Alpen Diamonds	R392F021FG	G392-314-021-10-F
60031839	G392-314-021-10.0-M, Alpen Diamonds	R392C021FG	G392-314-021-10-ML
60031841	G6051-314-018-06.0-M, Alpen Diamonds	R6051C018FG	G6051-314-018-6-ML
60031843	G6052-314-018-08.0-C, Alpen Diamonds	R6051SC018FG	G6051-314-018-8-MLX
60031845	G6052-314-018-08.0-M, Alpen Diamonds	R6052C018FG	G6052-314-018-8-ML
60031847	G6055-314-018-08.0-C, Alpen Diamonds	R6055SC018FG	G6055-314-018-8-MLX
60031849	G801-314-009-00.9-M, Alpen Diamonds	R801C009FG	G801-314-009-ML
60031851	G801-314-010-01.0-M, Alpen Diamonds	R801C010FG	G801-314-010-ML
60031853	G801-314-012-01.2-M, Alpen Diamonds	R801C012FG	G801-314-012-ML
60031857	G801-314-014-01.4-F, Alpen Diamonds	R801F014FG	G801-314-014-F
60031859	G801-314-014-01.4-M, Alpen Diamonds	R801C014FG	G801-314-014-ML
60031865	G801-314-016-01.6-M, Alpen Diamonds	R801C016FG	G801-314-016-ML
60031871	G801-314-018-01.8-F, Alpen Diamonds	R801F018FG	G801-314-018-F
60031873	G801-314-018-01.8-M, Alpen Diamonds	R801C018FG	G801-314-018-ML
60031877	G801-314-020-02.0-F, Alpen Diamonds	R801F021FG	G801-314-020-F
60031879	G801-314-020-02.0-M, Alpen Diamonds	R801C021FG	G801-314-020-ML
60031881	G801-314-023-02.3-C, Alpen Diamonds	R801SC023FG	G801-314-023-MLX
60031883	G801-314-023-02.3-EF, Alpen Diamonds	R801EF023FG	G801-314-023-XF
60031885	G801-314-023-02.3-F, Alpen Diamonds	R801F023FG	G801-314-023-F
60031887	G801-314-023-02.3-M, Alpen Diamonds	R801C023FG	G801-314-023-ML
60031891	G801-314-027-02.7-M, Alpen Diamonds	R801C027FG	G801-314-027-ML
60031893	G801-314-029-02.9-F, Alpen Diamonds	R801F029FG	G801-314-029-F
60031895	G801-314-029-02.9-M, Alpen Diamonds	R801C029FG	G801-314-029-ML
60031897	G801M-314-009-00.9-F, Alpen Diamonds	R801MF009FG	G801M-314-009-0.6-F
60031899	G801Z-314-012-01.2-M, Alpen Diamonds	ZREX801M012FG	ZR801-314-012-M
60031901	G801Z-314-014-01.4-M, Alpen Diamonds	ZREX801M014FG	ZR801-314-014-M
60031903	G801Z-314-018-01.8-M, Alpen Diamonds	ZREX801M018FG	ZR801-314-018-M
60031911	G805-314-012-00.9-M, Alpen Diamonds	R805C012FG	G805-314-012-0.9-ML
60031913	G805-314-014-01.4-M, Alpen Diamonds	R805C014FG	G805-314-014-1.4-ML
60031915	G805-314-016-01.4-M, Alpen Diamonds	R805C016FG	G805-314-016-1.4-ML
60031919	G806-314-012-02.2-M, Alpen Diamonds	R806C012FG	G806-314-012-2.2-ML
60031921	G806-314-014-02.2-M, Alpen Diamonds	R806C014FG	G806-314-014-2.2-ML
60031923	G806-314-016-02.5-M, Alpen Diamonds	R806C016FG	G806-314-016-2.5-ML
60031927	G807-314-012-03.5-M, Alpen Diamonds	R807C012FG	G807-314-012-3.5-ML
60031929	G807-314-016-04.0-M, Alpen Diamonds	R807C016FG	G807-314-016-4-ML
60031931	G807-314-018-05.0-M, Alpen Diamonds	R807C018FG	G807-314-018-5-ML
60031933	G811-314-031-04.2-M, Alpen Diamonds	R811C031FG	G811-314-031-4.2-ML
60031935	G811-314-033-05.0-C, Alpen Diamonds	R811SC033FG	G811-314-033-5-MLX
60031937	G811-314-033-05.0-M, Alpen Diamonds	R811C033FG	G811-314-033-5-ML
60031939	G811-314-037-07.0-C, Alpen Diamonds	R811SC037FG	G811-314-037-7-MLX
60031941	G811-314-037-07.0-M, Alpen Diamonds	R811C037FG	G811-314-037-7-ML
60031943	G811-314-048-08.0-M, Alpen Diamonds	R811C048FG	G811-314-048-8-ML
60031945	G811T-314-033-05.0-C, Alpen Diamonds	R811TC033FG	G811-314-033-5-TSX
60031949	G813-314-014-01.5-M, Alpen Diamonds	R813C014FG	G813-314-014-1.5-ML
60031951	G813-314-018-02.0-M, Alpen Diamonds	R813C018FG	G813-314-018-2-ML
60031955	G818-314-050-00.7-M, Alpen Diamonds	R818C050FG	G818-314-050-0.7-ML
60031961	G828-314-030-01.0-M, Alpen Diamonds	R828C030FG	G828-314-030-1.0-ML
60031963	G830-313-010-02.7-M, Alpen Diamonds	R830C0105S	G830-313-010-2.7-MLSS
60031967	G830-314-010-02.7-M, Alpen Diamonds	R830C010FG	G830-314-010-2.7-ML
60031969	G830-314-012-02.7-M, Alpen Diamonds	R830C012FG	G830-314-012-2.7-ML
60031973	G830-314-016-02.7-M, Alpen Diamonds	R830C016FG	G830-314-016-2.7-ML
60031981	G830L-314-012-04.0-M, Alpen Diamonds	R830LC012FG	G830L-314-012-4-ML
60031985	G830L-314-014-05.0-M, Alpen Diamonds	R830LC014FG	G830L-314-014-5-ML
60031987	G830L-314-016-05.0-C, Alpen Diamonds	R830LC016FG	G830L-314-016-5-ML
60031995	G834-314-016-06.8-M, Alpen Diamonds	R834C016FG	G834-314-016-6.8-ML
60031997	G834-314-021-06.8-M, Alpen Diamonds	R834C021FG	G834-314-021-6.8-ML
60031999	G835-313-010-04.0-M, Alpen Diamonds	R835C0105S	G835-313-010-4-MLSS
60032001	G835-314-008-03.0-M, Alpen Diamonds	R835C008FG	G835-314-008-3-ML
60032003	G835-314-010-04.0-M, Alpen Diamonds	R835C010FG	G835-314-010-4-ML
60032007	G835-314-012-04.0-M, Alpen Diamonds	R835C012FG	G835-314-012-4-ML
60032009	G835-314-014-04.0-M, Alpen Diamonds	R835C014FG	G835-314-014-4-ML
60032013	G835R-313-010-04.0-M, Alpen Diamonds	R835RC0105S	G835R-313-010-4-MLSS
60032015	G835R-314-008-03.0-M, Alpen Diamonds	R835RC008FG	G835R-314-008-3-ML
60032017	G835R-314-010-04.0-M, Alpen Diamonds	R835RC010FG	G835R-314-010-4-ML
60032021	G835R-314-012-04.0-M, Alpen Diamonds	R835RC012FG	G835R-314-012-4-ML
60032023	G835R-314-014-04.0-M, Alpen Diamonds	R835RC014FG	G835R-314-014-4-ML
60032025	G836-314-012-06.0-F, Alpen Diamonds	R836F012FG	G836-314-012-6-ML
60032027	G836-314-012-06.0-M, Alpen Diamonds	R836C012FG	G836-314-012-6-ML
60032029	G836-314-014-06.0-C, Alpen Diamonds	R836SC014FG	G836-314-014-6-MLX
60032031	G836-314-014-06.0-M, Alpen Diamonds	R836C014FG	G836-314-014-6-ML
60032037	G836R-314-010-06.0-M, Alpen Diamonds	R836RC010FG	G836R-314-010-6-ML
60032039	G836R-314-012-06.0-M, Alpen Diamonds	R836RC012FG	G836R-314-012-6-ML

REF new	Description new	REF old	Description old
60032041	G836R-314-014-06.0-F, Alpen Diamonds	R836RF014FG	G836R-314-014-6-F
60032043	G836R-314-014-06.0-M, Alpen Diamonds	R836RC014FG	G836R-314-014-6-ML
60032045	G836T-314-014-06.0-M, Alpen Diamonds	R836TC014FG	G836-314-014-6-TSX
60032047	G836T-314-016-06.0-C, Alpen Diamonds	R836TC016FG	G836-314-016-6-TSX
60032049	G837-314-012-07.0-M, Alpen Diamonds	R837C012FG	G837-314-012-7-ML
60032051	G837-314-012-08.0-M, Alpen Diamonds	R837RC012FG	G837-314-012-8-ML
60032053	G837-314-014-08.0-M, Alpen Diamonds	R837C014FG	G837-314-014-8-ML
60032055	G837-314-016-08.0-C, Alpen Diamonds	R837SC016FG	G837-314-016-8-MLX
60032057	G837-314-016-08.0-M, Alpen Diamonds	R837C016FG	G837-314-016-8-ML
60032061	G837R-314-012-08.0-M, Alpen Diamonds	R837RC012FG	G837R-314-012-8-ML
60032063	G837R-314-014-08.0-F, Alpen Diamonds	R837RF014FG	G837R-314-014-8-F
60032065	G837R-314-014-08.0-M, Alpen Diamonds	R837RC014FG	G837R-314-014-8-ML
60032067	G837R-314-016-08.0-M, Alpen Diamonds	R837RC016FG	G837R-314-016-8-ML
60032069	G837T-314-014-08.0-M, Alpen Diamonds	R837TC014FG	G837-314-014-8-TSX
60032071	G837T-314-016-08.0-C, Alpen Diamonds	R837TC016FG	G837-314-016-8-TSX
60032075	G838-314-010-04.0-M, Alpen Diamonds	R838C010FG	G838-314-010-4-ML
60032077	G838-314-012-04.0-M, Alpen Diamonds	R838C012FG	G838-314-012-4-ML
60032081	G838M-314-007-02.5-F, Alpen Diamonds	R838MF007FG	G838M-314-007-2.4-F
60032083	G839-314-010-00.2-M, Alpen Diamonds	R839C010FG	G839-314-010-ML
60032085	G839-314-012-00.2-F, Alpen Diamonds	R839F012FG	G839-314-012-F
60032087	G839-314-012-00.2-M, Alpen Diamonds	R839C012FG	G839-314-012-ML
60032089	G839-314-014-00.2-F, Alpen Diamonds	R839F014FG	G839-314-014-F
60032091	G839-314-014-00.2-M, Alpen Diamonds	R839C014FG	G839-314-014-ML
60032095	G845-314-010-04.0-M, Alpen Diamonds	R845C010FG	G845-314-010-4-ML
60032097	G845-314-012-04.0-M, Alpen Diamonds	R845C012FG	G845-314-012-4-ML
60032103	G845-314-016-04.0-M, Alpen Diamonds	R845C016FG	G845-314-016-4-ML
60032105	G845R-314-016-04.0-M, Alpen Diamonds	R845RC016FG	G845R-314-016-4-ML
60032107	G845R-314-025-04.0-F, Alpen Diamonds	R845RF025FG	G845R-314-025-4-F
60032109	G845R-314-025-04.0-M, Alpen Diamonds	R845RC025FG	G845R-314-025-4-ML
60032111	G846-314-016-06.0-M, Alpen Diamonds	R846C016FG	G846-314-016-6-ML
60032113	G846-314-016-07.0-M, Alpen Diamonds	R846F016FG	G846-314-016-7-ML
60032115	G846-314-018-07.0-M, Alpen Diamonds	R846C018FG	G846-314-018-7-ML
60032117	G846-314-023-07.0-M, Alpen Diamonds	R846C023FG	G846-314-023-7-ML
60032119	G846R-314-016-06.0-F, Alpen Diamonds	R846RF016FG	G846R-314-016-6-F
60032121	G846R-314-016-06.0-M, Alpen Diamonds	R846RC016FG	G846R-314-016-6-ML
60032123	G847-313-014-08.0-M, Alpen Diamonds	R847SC0145S	G847-313-014-8-MLXSS
60032125	G847-314-012-08.0-M, Alpen Diamonds	R847C012FG	G847-314-012-8-ML
60032127	G847-314-014-08.0-F, Alpen Diamonds	R847F014FG	G847-314-014-8-F
60032129	G847-314-014-08.0-M, Alpen Diamonds	R847C014FG	G847-314-014-8-ML
60032131	G847-314-016-08.0-C, Alpen Diamonds	R847SC016FG	G847-314-016-8-MLX
60032133	G847-314-016-08.0-F, Alpen Diamonds	R847F016FG	G847-314-016-8-F
60032135	G847-314-016-08.0-M, Alpen Diamonds	R847C016FG	G847-314-016-8-ML
60032137	G847-314-018-08.0-C, Alpen Diamonds	R847SC018FG	G847-314-018-8-MLX
60032139	G847-314-018-08.0-F, Alpen Diamonds	R847F018FG	G847-314-018-8-F
60032141	G847-314-018-08.0-M, Alpen Diamonds	R847C018FG	G847-314-018-8-ML
60032143	G847-314-025-08.0-C, Alpen Diamonds	R847SC025FG	G847-314-025-8-MLX
60032145	G847-314-025-08.0-M, Alpen Diamonds	R847C025FG	G847-314-025-8-ML
60032147	G847R-314-016-08.0-C, Alpen Diamonds	R847RC016FG	G847R-314-016-8-ML
60032149	G847R-314-016-08.0-F, Alpen Diamonds	R847RF016FG	G847R-314-016-8-F
60032155	G847T-314-016-08.0-C, Alpen Diamonds	R847TC016FG	G847-314-016-8-TSX
60032157	G847T-314-018-08.0-C, Alpen Diamonds	R847TC018FG	G847-314-018-8-TSX
60032159	G848-314-012-10.0-M, Alpen Diamonds	R848C012FG	G848-314-012-10-ML
60032163	G848-314-014-10.0-M, Alpen Diamonds	R848C014FG	G848-314-014-10-ML
60032165	G848-314-016-10.0-C, Alpen Diamonds	R848C016FG	G848-314-016-10-ML
60032171	G848-314-018-10.0-C, Alpen Diamonds	R848SC018FG	G848-314-018-10-MLX
60032173	G848-314-018-10.0-M, Alpen Diamonds	R848C018FG	G848-314-018-10-ML
60032175	G848-314-023-10.0-C, Alpen Diamonds	R848SC023FG	G848-314-023-10-MLX
60032177	G848-314-023-10.0-F, Alpen Diamonds	R848F023FG	G848-314-023-10-F
60032179	G848-314-023-10.0-M, Alpen Diamonds	R848C023FG	G848-314-023-10-ML
60032181	G848R-314-016-10.0-F, Alpen Diamonds	R850RF018FG	G848R-314-016-10-F
60032183	G848R-314-018-10.0-M, Alpen Diamonds	R850RC018FG	G848R-314-018-10-ML
60032185	G848T-314-016-10.0-C, Alpen Diamonds	R848TC016FG	G848-314-016-10-TSX
60032187	G848T-314-018-10.0-C, Alpen Diamonds	R848TC018FG	G848-314-018-10-TSX
60032189	G849-314-016-04.0-M, Alpen Diamonds	R849C016FG	G849-314-016-4-ML
60032191	G849-314-025-04.0-C, Alpen Diamonds	R849SC025FG	G849-314-025-4-MLX
60032193	G849-314-025-04.0-M, Alpen Diamonds	R849C025FG	G849-314-025-4-ML
60032195	G850-314-012-10.0-F, Alpen Diamonds	R850F012FG	G850-314-012-10-F
60032197	G850-314-012-10.0-M, Alpen Diamonds	R850C012FG	G850-314-012-10-ML
60032199	G850-314-014-10.0-C, Alpen Diamonds	R850SC014FG	G850-314-014-10-MLX
60032201	G850-314-014-10.0-F, Alpen Diamonds	R850F014FG	G850-314-014-10-F
60032203	G850-314-014-10.0-M, Alpen Diamonds	R850C014FG	G850-314-014-10-ML
60032205	G850-314-016-10.0-C, Alpen Diamonds	R850SC016FG	G850-314-016-10-MLX
60032207	G850-314-016-10.0-F, Alpen Diamonds	R850F016FG	G850-314-016-10-F
60032209	G850-314-016-10.0-M, Alpen Diamonds	R850C016FG	G850-314-016-10-ML
60032211	G850-314-018-08.0-C, Alpen Diamonds	R850SC018FG	G850-314-018-8-MLX
60032213	G850-314-018-08.0-M, Alpen Diamonds	R850RC018FG	G850-314-018-8-ML
60032217	G850-314-018-10.0-C, Alpen Diamonds	R850T018C018FG	G850-314-018-10-MLX
60032219	G850-314-018-10.0-F, Alpen Diamonds	R850F018FG	G850-314-018-10-F

REF new	Description new	REF old	Description old
60032221	G850-314-018-10.0-M, Alpen Diamonds	R85010C018FG	G850-314-018-10-ML
60032223	G850-314-023-10.0-M, Alpen Diamonds	R850C023FG	G850-314-023-10-ML
60032225	G850L-314-016-12.0-F, Alpen Diamonds	R850L016FG	G850L-314-016-12-F
60032227	G850L-314-016-12.0-M, Alpen Diamonds	R850LC016FG	G850L-314-016-12-ML
60032229	G850L-314-018-12.0-M, Alpen Diamonds	R850LC018FG	G850L-314-018-12-ML
60032233	G850T-314-014-10.0-M, Alpen Diamonds	R850T0C14FG	G850-314-014-10-TSX
60032235	G850T-314-016-10.0-C, Alpen Diamonds	R850T0C16FG	G850-314-016-10-TSX
60032237	G850T-314-018-08.0-C, Alpen Diamonds	R850T0C18FG	G850-314-018-8-TSX
60032239	G850T-314-018-10.0-C, Alpen Diamonds	R85010T0C18FG	G850-314-018-10-TSX
60032241	G850Z-314-016-10.0-M, Alpen Diamonds	ZREX850M016FG	ZR850-314-016-10-M
60032243	G850Z-314-018-10.0-M, Alpen Diamonds	ZREX850M018FG	ZR850-314-018-10-M
60032245	G851-314-012-08.0-M, Alpen Diamonds	R851C012FG	G851-314-012-8-ML
60032247	G851-314-016-08.0-M, Alpen Diamonds	R851C016FG	G851-314-016-8-ML
60032249	G852-314-012-06.0-F, Alpen Diamonds	R852E0F12FG	G852-314-012-6-XF
60032251	G852-314-012-06.0-F, Alpen Diamonds	R852F012FG	G852-314-012-6-F
60032253	G852-314-012-06.0-M, Alpen Diamonds	R852C012FG	G852-314-012-6-ML
60032255	G852-314-016-06.0-F, Alpen Diamonds	R852E0F16FG	G852-314-016-6-XF
60032259	G852-314-016-06.0-M, Alpen Diamonds	R852C016FG	G852-314-016-6-ML
60032261	G853-314-008-03.5-F, Alpen Diamonds	R853E008FG	G853-314-008-3.5-XF
60032263	G853-314-008-03.5-M, Alpen Diamonds	R853C008FG	G853-314-008-3.5-ML
60032265	G855-313-014-06.0-M, Alpen Diamonds	R855C014SS	G855-313-014-6-MLSS
60032267	G855-314-010-06.0-M, Alpen Diamonds	R855C010FG	G855-314-010-6-ML
60032269	G855-314-012-06.0-M, Alpen Diamonds	R855C012FG	G855-314-012-6-ML
60032271	G855-314-014-06.0-M, Alpen Diamonds	R855C014FG	G855-314-014-6-ML
60032273	G856-313-016-07.0-M, Alpen Diamonds	R856C016SS	G856-313-016-7-MLSS
60032275	G856-313-018-08.0-M, Alpen Diamonds	R856C018SS	G856-313-018-8-MLSS
60032277	G856-314-012-08.0-F, Alpen Diamonds	R856F012FG	G856-314-012-8-F
60032279	G856-314-012-08.0-M, Alpen Diamonds	R856C012FG	G856-314-012-8-ML
60032281	G856-314-014-08.0-F, Alpen Diamonds	R856F014FG	G856-314-014-8-F
60032283	G856-314-014-08.0-M, Alpen Diamonds	R856C014FG	G856-314-014-8-ML
60032285	G856-314-014-09.0-M, Alpen Diamonds	R856L9C014FG	G856-314-014-9-ML
60032287	G856-314-016-07.0-C, Alpen Diamonds	R856C016FG	G856-314-016-7-ML
60032289	G856-314-016-07.0-F, Alpen Diamonds	R856F016FG	G856-314-016-7-F
60032293	G856-314-016-08.0-C, Alpen Diamonds	R856S0C016FG	G856-314-016-8-MLX
60032295	G856-314-016-08.0-F, Alpen Diamonds	R856S0F016FG	G856-314-016-8-F
60032297	G856-314-016-08.0-M, Alpen Diamonds	R856S0C016FG	G856-314-016-8-ML
60032305	G856-314-018-08.0-C, Alpen Diamonds	R856S0C18FG	G856-314-018-8-MLX
60032307	G856-314-018-08.0-F, Alpen Diamonds	R856F018FG	G856-314-018-8-F
60032309	G856-314-018-08.0-M, Alpen Diamonds	R856C018FG	G856-314-018-8-ML
60032311	G856-314-018-09.0-C, Alpen Diamonds	R856LSC018FG	G856-314-018-9-MLX
60032319	G856-314-021-08.0-C, Alpen Diamonds	R856S0C21FG	G856-314-021-8-MLX
60032321	G856-314-021-08.0-M, Alpen Diamonds	R856C021FG	G856-314-021-8-ML
60032323	G856-314-021-09.0-C, Alpen Diamonds	R856LSC021FG	G856-314-021-9-MLX
60032325	G856-314-021-09.0-F, Alpen Diamonds	R856LF021FG	G856-314-021-9-F
60032327	G856-314-021-09.0-M, Alpen Diamonds	R856LC021FG	G856-314-021-9-ML
60032329	G856-314-025-07.0-C, Alpen Diamonds	R855SC025FG	G856-314-025-7-MLX
60032333	G856-314-025-07.0-M, Alpen Diamonds	R855C025FG	G856-314-025-7-ML
60032335	G856-314-025-08.0-C, Alpen Diamonds	R856S0C25FG	G856-314-025-8-MLX
60032337	G856-314-025-08.0-M, Alpen Diamonds	R856C025FG	G856-314-025-8-ML
60032341	G856L-314-014-10.0-M, Alpen Diamonds	R856L10C014FG	G856L-314-014-10-ML
60032343	G856L-314-016-10.0-M, Alpen Diamonds	R856LC016FG	G856-314-016-9-MLX
60032345	G856L-314-016-10.0-M, Alpen Diamonds	R856C016FG	G856-314-016-9-ML
60032349	G856L-314-018-10.0-M, Alpen Diamonds	R856LC018FG	G856L-314-018-9-ML
60032357	G856LT-314-016-10.0-C, Alpen Diamonds	R856LTC016FG	G856L-314-016-10-TSX
60032359	G856LT-314-018-11.0-C, Alpen Diamonds	R856LTC018FG	G856L-314-018-11-TSX
60032361	G856N-314-014-08.0-M, Alpen Diamonds	R856NC014FG	G856N-314-014-8-ML
60032367	G856T-314-014-08.0-M, Alpen Diamonds	R856T0C14FG	G856-314-014-8-TSX
60032369	G856T-314-018-08.0-C, Alpen Diamonds	R856T0C18FG	G856-314-018-8-TSX
60032371	G856T-314-021-09.0-C, Alpen Diamonds	R856T0C21FG	G856-314-021-9-TSX
60032373	G856T-314-025-07.0-C, Alpen Diamonds	R856T0C25FG	G856-314-025-7-TSX
60032379	G856Z-314-018-08.0-M, Alpen Diamonds	ZREX856M018FG	ZR856-314-018-8-M
60032381	G858-313-012-08.0-M, Alpen Diamonds	R858C012SS	G858-313-012-8-MLSS
60032383	G858-314-012-08.0-EF, Alpen Diamonds	R858E0F12FG	G858-314-012-8-XF
60032385	G858-314-012-08.0-F, Alpen Diamonds	R858F012FG	G858-314-012-8-F
60032387	G858-314-012-08.0-M, Alpen Diamonds	R858C012FG	G858-314-012-8-ML
60032389	G858-314-014-08.0-EF, Alpen Diamonds	R858E0F14FG	G858-314-014-8-XF
60032391	G858-314-014-08.0-F, Alpen Diamonds	R858F014FG	G858-314-014-8-F
60032393	G858-314-014-08.0-M, Alpen Diamonds	R858C014FG	G858-314-014-8-ML
60032395	G858-314-016-08.0-M, Alpen Diamonds	R858C016FG	G858-314-016-8-ML
60032397	G859-314-010-10.0-EF, Alpen Diamonds	R859E0F10FG	G859-314-010-10-XF
60032399	G859-314-010-10.0-F, Alpen Diamonds	R859F010FG	G859-314-010-10-F
60032403	G859-314-010-10.0-M, Alpen Diamonds	R859C010FG	G859-314-010-10-ML
60032409	G859-314-014-10.0-EF, Alpen Diamonds	R859E0F14FG	G859-314-014-10-XF
60032411	G859-314-014-10.0-F, Alpen Diamonds	R859F014FG	G859-314-014-10-F
60032415	G859-314-014-10.0-M, Alpen Diamonds	R859C014FG	G859-314-014-10-ML
60032419	G859-314-018-10.0-EF, Alpen Diamonds	R859E0F18FG	G859-314-018-10-XF
60032421	G859-314-018-10.0-F, Alpen Diamonds	R859F018FG	G859-314-018-10-F
60032423	G859-314-018-10.0-M, Alpen Diamonds	R859C018FG	G859-314-018-10-ML

REF new	Description new	REF old	Description old
60032425	G860-314-010-04.0-M, Alpen Diamonds	R860C010FG	G860-314-010-4-ML
60032427	G860-314-012-05.0-EF, Alpen Diamonds	R8605EF012FG	G860-314-012-5-XF
60032429	G860-314-012-05.0-F, Alpen Diamonds	R8605F012FG	G860-314-012-5-F
60032431	G860-314-012-05.0-M, Alpen Diamonds	R8605C012FG	G860-314-012-5-ML
60032433	G860-314-012-06.5-F, Alpen Diamonds	R86065F012FG	G860-314-012-6.5-F
60032435	G860-314-012-06.5-M, Alpen Diamonds	R86065C012FG	G860-314-012-6.5-ML
60032439	G860-314-014-05.0-F, Alpen Diamonds	R860C014FG	G860-314-014-5-F
60032441	G860-314-014-05.0-M, Alpen Diamonds	R860C014FG	G860-314-014-5-ML
60032443	G862-313-010-08.0-M, Alpen Diamonds	R862C010SS	G862-313-010-8-MLSS
60032445	G862-314-010-08.0-EF, Alpen Diamonds	R862E0F10FG	G862-314-010-8-XF
60032447	G862-314-010-08.0-F, Alpen Diamonds	R862F010FG	G862-314-010-8-F
60032449	G862-314-010-08.0-M, Alpen Diamonds	R862C010FG	G862-314-010-8-ML
60032451	G862-314-012-08.0-EF, Alpen Diamonds	R862E0F12FG	G862-314-012-8-XF
60032453	G862-314-012-08.0-F, Alpen Diamonds	R862F012FG	G862-314-012-8-F
60032455	G862-314-012-08.0-M, Alpen Diamonds	R862C012FG	G862-314-012-8-ML
60032457	G862-314-014-08.0-F, Alpen Diamonds	R862F014FG	G862-314-014-8-F
60032459	G862-314-014-08.0-M, Alpen Diamonds	R862C014FG	G862-314-014-8-ML
60032461	G862-314-016-08.0-EF, Alpen Diamonds	R862E0F16FG	G862-314-016-8-XF
60032463	G862-314-016-08.0-F, Alpen Diamonds	R862F016FG	G862-314-016-8-F
60032465	G862-314-016-08.0-M, Alpen Diamonds	R862C016FG	G862-314-016-8-ML
60032467	G862Z-314-014-08.0-M, Alpen Diamonds	ZREX86ZM014FG	ZR862-314-014-8-M
60032469	G863-314-012-10.0-EF, Alpen Diamonds	R863E0F12FG	G863-314-012-10-XF
60032471	G863-314-012-10.0-F, Alpen Diamonds	R863F012FG	G863-314-012-10-F
60032473	G863-314-012-10.0-M, Alpen Diamonds	R863C012FG	G863-314-012-10-ML
60032475	G863-314-016-10.0-EF, Alpen Diamonds	R863E0F16FG	G863-314-016-10-XF
60032477	G863-314-016-10.0-F, Alpen Diamonds	R863F016FG	G863-314-016-10-F
60032479	G863-314-016-10.0-M, Alpen Diamonds	R863C016FG	G863-314-016-10-ML
60032481	G863Z-314-014-10.0-M, Alpen Diamonds	ZREX86ZM014FG	ZR863-314-014-10-M
60032483	G864-314-014-12.0-F, Alpen Diamonds	R864F014FG	G864-314-014-12-F
60032485	G864-314-014-12.0-M, Alpen Diamonds	R864C014FG	G864-314-014-12-ML
60032487	G877-314-010-06.0-M, Alpen Diamonds	R877C010FG	G877-314-010-6-ML
60032489	G877-314-012-06.0-M, Alpen Diamonds	R877C012FG	G877-314-012-6-ML
60032491	G877K-314-012-06.0-M, Alpen Diamonds	R877KC012FG	G877K-314-012-6-ML
60032495	G877K-314-014-06.0-M, Alpen Diamonds	R877KC014FG	G877K-314-014-6-ML
60032497	G877K-314-016-06.0-M, Alpen Diamonds	R877KC016FG	G877K-314-016-6-ML
60032499	G877K-314-018-06.0-M, Alpen Diamonds	R877KC018FG	G877K-314-018-6-ML
60032501	G878-314-010-08.0-F, Alpen Diamonds	R878F010FG	G878-314-010-8-F
60032503	G878-314-010-08.0-M, Alpen Diamonds	R878C010FG	G878-314-010-8-ML
60032505	G878-314-012-08.0-F, Alpen Diamonds	R878F012FG	G878-314-012-8-F
60032507	G878-314-012-08.0-M, Alpen Diamonds	R878C012FG	G878-314-012-8-ML
60032509	G878-314-014-08.0-M, Alpen Diamonds	R878C014FG	G878-314-014-8-ML
60032511	G878-314-016-08.0-M, Alpen Diamonds	R878C016FG	G878-314-016-8-ML
60032513	G878K-314-010-08.0-F, Alpen Diamonds	R878KF010FG	G878K-314-010-8-F
60032515	G878K-314-010-08.0-M, Alpen Diamonds	R878KC010FG	G878K-314-010-8-ML
60032517	G878K-314-012-08.0-F, Alpen Diamonds	R878KF012FG	G878K-314-012-8-F
60032519	G878K-314-012-08.0-M, Alpen Diamonds	R878KC012FG	G878K-314-012-8-ML
60032521	G878K-314-014-08.0-F, Alpen Diamonds	R878KF014FG	G878K-314-014-8-F
60032523	G878K-314-014-08.0-M, Alpen Diamonds	R878KC014FG	G878K-314-014-8-ML
60032525	G878K-314-016-08.0-C, Alpen Diamonds	R878KSC016FG	G878K-314-016-8-MLX
60032529	G878K-314-016-08.0-M, Alpen Diamonds	R878KC016FG	G878K-314-016-8-ML
60032531	G878K-314-018-08.0-C, Alpen Diamonds	R878KSC018FG	G878K-314-018-8-MLX
60032533	G878K-314-018-08.0-F, Alpen Diamonds	R878KF018FG	G878K-314-018-8-F
60032535	G878K-314-018-08.0-M, Alpen Diamonds	R878KC018FG	G878K-314-018-8-ML
60032537	G878K-314-021-08.0-C, Alpen Diamonds	R878KSC021FG	G878K-314-021-8-MLX
60032539	G878K-314-021-08.0-M, Alpen Diamonds	R878KC021FG	G878K-314-021-8-ML
60032549	G878KT-314-016-08.0-C, Alpen Diamonds	R878KTC016FG	G878K-314-016-8-TSX
60032553	G879-314-012-10.0-M, Alpen Diamonds	R879C012FG	G879-314-012-10-ML
60032555	G879-314-014-10.0-C, Alpen Diamonds	R879C014FG	G879-314-014-10-ML
60032561	G879-314-016-10.0-M, Alpen Diamonds	R879C016FG	G879-314-016-10-ML
60032563	G879K-314-012-10.0-M, Alpen Diamonds	R879KC012FG	G879K-314-012-10-ML
60032565	G879K-314-014-10.0-M, Alpen Diamonds	R879KC014FG	G879K-314-014-10-ML
60032567	G879K-314-016-10.0-C, Alpen Diamonds	R879KSC016FG	G879K-314-016-10-MLX
60032569	G879K-314-016-10.0-M, Alpen Diamonds	R879KC016FG	G879K-314-016-10-ML
60032571	G879K-314-018-10.0-C, Alpen Diamonds	R879KSC018FG	G879K-314-018-10-MLX
60032575	G879K-314-018-10.0-M, Alpen Diamonds	R879KC018FG	G879K-314-018-10-ML
60032577	G879K-314-021-10.0-C, Alpen Diamonds	R879KSC021FG	G879K-314-021-10-MLX
60032579	G879K-314-021-10.0-M, Alpen Diamonds	R879KC021FG	G879K-314-021-10-ML
60032581	G879K-314-023-10.0-M, Alpen Diamonds	R879KC023FG	G879K-314-023-10-ML
60032583	G879KT-314-014-10.0-M, Alpen Diamonds	R879KTC014FG	G879K-314-014-10-TSX
60032585	G879KT-314-018-10.0-C, Alpen Diamonds	R879KTC018FG	G879K-314-018-10-TSX
60032589	G880-314-012-06.0-M, Alpen Diamonds	R880C012FG	G880-314-012-6-ML
60032591	G880-314-014-06.0-M, Alpen Diamonds	R880C014FG	G880-314-014-6-ML
60032595	G881-314-010-08.0-M, Alpen Diamonds	R881C010FG	G881-314-010-8-ML
60032601	G881-314-012-08.0-M, Alpen Diamonds	R881C012FG	G881-314-012-8-ML
60032603	G881-314-014-08.0-C, Alpen Diamonds	R881SC014FG	G881-314-014-8-MLX
60032607	G881-314-014-08.0-M, Alpen Diamonds	R881C014FG	G881-314-014-8-ML
60032609	G881-314-016-08.0-C, Alpen Diamonds	R881SC016FG	G881-314-016-8-MLX
60032611	G881-314-016-08.0-F, Alpen Diamonds	R881F016FG	G881-314-016-8-F
60032613	G881-314-016-08.0-M, Alpen Diamonds	R881C016FG	G881-314-016-8-ML

REF new	Description new	REF old	Description old
60032615	G881T-314-014-08.0-M, Alpen Diamonds	R881TC014FG	G881-314-014-8-TSX
60032617	G881TC-314-016-08.0-C, Alpen Diamonds	R881TC014FG	G881-314-016-8-TSX
60032619	G881Z-314-014-08.0-M, Alpen Diamonds	ZREX881M014FG	ZR881-314-014-8-M
60032621	G881Z-314-016-08.0-C, Alpen Diamonds	ZREX881C016FG	ZR881-314-016-8-C
60032627	G882-314-014-10.0-M, Alpen Diamonds	R882C014FG	G882-314-014-10-ML
60032629	G882T-314-014-10.0-M, Alpen Diamonds	R882TC014FG	G882-314-014-10-TSX
60032631	G884-314-012-06.0-M, Alpen Diamonds	R884C012FG	G884-314-012-6-ML
60032633	G885-313-012-08.0-M, Alpen Diamonds	R885C012SS	G885-313-012-8-MLS5
60032635	G885-314-012-08.0-F, Alpen Diamonds	R885F012FG	G885-314-012-8-F
60032637	G885-314-012-08.0-M, Alpen Diamonds	R885C012FG	G885-314-012-8-ML
60032639	G885-314-014-08.0-C, Alpen Diamonds	R885C014FG	G885-314-014-8-MLX
60032641	G885-314-014-08.0-M, Alpen Diamonds	R885C014FG	G885-314-014-8-ML
60032643	G886-314-012-10.0-F, Alpen Diamonds	R886F012FG	G886-314-012-10-F
60032645	G886-314-012-10.0-M, Alpen Diamonds	R886C012FG	G886-314-012-10-ML
60032651	G886-314-016-10.0-C, Alpen Diamonds	R886C014FG	G886-314-016-10-MLX
60032653	G886-314-016-10.0-M, Alpen Diamonds	R886C016FG	G886-314-016-10-ML
60032655	G889-314-009-03.5-EF, Alpen Diamonds	R889EF009FG	G889-314-009-3.5-XF
60032657	G889-314-009-03.5-F, Alpen Diamonds	R889F009FG	G889-314-009-3.5-F
60032659	G889-314-009-03.5-M, Alpen Diamonds	R889C009FG	G889-314-009-3.5-ML
60032661	G889M-314-007-02.5-F, Alpen Diamonds	R889M007FG	G889M-314-007-2.4-F
60032663	G899-314-027-07.0-M, Alpen Diamonds	R899C027FG	G899-314-027-7-ML
60032665	G905-314-023-02.7-F, Alpen Diamonds	R905F023FG	G905-314-023-2.7-F
60032667	G905-314-027-02.9-F, Alpen Diamonds	R905F027FG	G905-314-027-2.9-F
60032669	G905-314-031-03.1-F, Alpen Diamonds	R905C031FG	G905-314-031-3.1-F
60032671	G909-314-040-01.0-M, Alpen Diamonds	R909I040FG	G909-314-040-1-ML
60032673	G909-314-040-01.5-C, Alpen Diamonds	R909I15SC040FG	G909-314-040-1.5-MLX
60032675	G909-314-040-01.5-M, Alpen Diamonds	R909I15C040FG	G909-314-040-1.5-ML
60032677	G909-314-055-02.0-M, Alpen Diamonds	R909C055FG	G909-314-055-2-ML
60032679	GET3-314-008-03.0-F, Alpen Diamonds	R132F008FG	G/ET3-314-008-3-F
60032681	GET4-314-010-04.0-F, Alpen Diamonds	R133F01FG	G/ET4-314-010-4-F
60032683	GET4-314-010-04.2-EF, Alpen Diamonds	R133EF01FG	G/ET4-314-010-4.2-XF
60032685	GET6-314-014-06.0-EF, Alpen Diamonds	R134EF014FG	G/ET6-314-014-6-XF
60032687	GET6-314-014-06.0-F, Alpen Diamonds	R134F014FG	G/ET6-314-014-6-F
60032689	GET9-314-014-09.0-EF, Alpen Diamonds	R135EF014FG	G/ET9-314-014-9-XF
60032691	GET9-314-014-09.0-F, Alpen Diamonds	R135F014FG	G/ET9-314-014-9-F
60032695	KK50-314-010-07.0-M, Alpen Diamonds	RK50C010FG	KS0-010-7-ML
60032697	KK51-314-012-07.0-M, Alpen Diamonds	RK51SC012FG	KS1-012-7-MLX
60032699	KK52-314-014-08.0-C, Alpen Diamonds	RK52SC014FG	KS2-014-8-MLX
60032701	KK54-314-024-04.2-C, Alpen Diamonds	RK54SC024FG	KS4-024-4.2-MLX
60032703	KK54-314-024-04.2-F, Alpen Diamonds	RK54F024FG	KS4-024-4.2-ML
60032705	KK55-314-031-06.2-C, Alpen Diamonds	RK55SC031FG	KS5-031-6.2-MLX
60032707	KK56-314-025-09.0-M, Alpen Diamonds	RK56C025FG	KS6-025-9.0-ML
60032709	KK57-314-054-01.0-M, Alpen Diamonds	RK57SC054FG	KS7-054-1-MLX