



Diamonds

ORDER INFO:
Please add prefix
70-D
EX: 70-D801009M

ROUND						
Code	801009	801012	801014	801016	801018	801022
Diameter (mm)	0.9	1.2	1.4	1.6	1.8	2.2
Head Length (mm)	0.9	1.2	1.4	1.6	1.8	2.2
Super Fine						
Fine						
Medium	M	M	M	M	M	M
Coarse		C	C	C	C	C
Super Coarse						

BARREL		
Code	811033	811037
Diameter (mm)	3.3	3.7
Head Length (mm)	5.5	5.5
Super Fine		
Fine		
Medium		
Coarse	C	C
Super Coarse		SC

INVERTED CONE						
Code	805012	805016	805018	807014	807016	807018
Diameter (mm)	1.2	1.6	1.8	1.4	1.6	1.8
Head Length (mm)	1.3	1.4	1.6	3.0	4.0	5.0
Super Fine						
Fine				F		
Medium				M	M	M
Coarse	C	C	C			
Super Coarse						

DOUBLE INVERTED CONE			
Code	806012	806014	813018
Diameter (mm)	1.2	1.4	1.8
Head Length (mm)	3.0	2.0	2.1
Super Fine			
Fine			
Medium		M	M
Coarse	C		
Super Coarse			



PEAR	
Code	830L012
Diameter (mm)	1.2
Head Length (mm)	3.9
Super Fine	
Fine	
Medium	M
Coarse	
Super Coarse	

WHEEL	
Code	818030
Diameter (mm)	3.0
Head Length (mm)	0.8
Super Fine	
Fine	
Medium	M
Coarse	
Super Coarse	

ROUND WHEEL		
Code	909037	909040
Diameter (mm)	3.7	4.0
Head Length (mm)	1.0	1.3
Super Fine		
Fine		
Medium	M	
Coarse	C	C
Super Coarse		

INTER PROXIMAL	
Code	392016
Diameter (mm)	1.6
Head Length (mm)	5.0
Super Fine	SF
Fine	F
Medium	
Coarse	C
Super Coarse	

GROSS REDUCTION				
Code	651016	652018	654016	655018
Diameter (mm)	1.6	1.8	1.6	1.8
Head Length (mm)	6.0	8.0	6.0	8.0
Super Fine				
Fine				
Medium				
Coarse	C	C	C	C
Super Coarse				

EGG/FOOTBALL							
Code	368016	368018	368023	369023	379016	379018	379023
Diameter (mm)	1.6	1.8	2.3	2.3	1.6	1.8	2.3
Head Length (mm)	3.3	4.5	5.3	5.1	3.4	3.6	4.4
Super Fine		SF	SF				SF
Fine	F	F	F	F	F	F	F
Medium		M	M			M	M
Coarse		C	C	C		C	C
Super Coarse			SC				

FLAT END CYLINDER											
Code	835008	835010	835012	835014	836012	836014	837010	837012	837014	837016	837018
Diameter (mm)	0.8	1.0	1.2	1.4	1.2	1.4	1.0	1.2	1.4	1.6	1.8
Head Length (mm)	4.0	4.0	4.0	4.0	5.6	6.0	8.0	8.0	8.0	8.0	8.0
Super Fine											
Fine											
Medium	M	M						M	M	M	
Coarse		C	C	C	C	C	C	C	C	C	C
Super Coarse											

TORPEDO CYLINDER		
Code	878012	878014
Diameter (mm)	1.2	1.4
Head Length (mm)	8.0	8.0
Super Fine		
Fine		
Medium	M	
Coarse	C	C
Super Coarse		

CYLINDER ROUNDED SHOULDER			
Code	836KR014	837KR012	837KR016
Diameter (mm)	1.4	1.2	1.6
Head Length (mm)	6.0	8.0	8.0
Super Fine			
Fine			
Medium			
Coarse	C	C	C
Super Coarse			

ROUND END CYLINDER		
Code	881014	881016
Diameter (mm)	1.4	1.6
Head Length (mm)	8.0	8.0
Super Fine		
Fine		
Medium		M
Coarse	C	C
Super Coarse		

SF: Super Fine
F: Fine
M: Medium
C: Coarse
SC: Super Coarse



Diamonds

FLAT END TAPER									
Code	846012	846018	847014	847016	847018	847019	848012	848016	848018
Diameter (mm)	1.2	1.8	1.4	1.6	1.8	1.9	1.2	1.6	1.8
Head Length (mm)	6.0	7.0	8.0	8.0	8.0	8.0	10.0	10.0	10.0
Super Fine									
Fine									
Medium				Ⓜ		Ⓜ		Ⓜ	
Coarse	C	C	C	C	C	C	C	C	C
Super Coarse									

POINTED CYLINDER				
Code	885010	885012	886012	886016
Diameter (mm)	1.0	1.2	1.2	1.6
Head Length (mm)	8.0	8.0	10.0	10.0
Super Fine				
Fine				
Medium		Ⓜ	Ⓜ	
Coarse	C	C	C	C
Super Coarse				

GINGIVAL CURETTAGE				
Code	878K014	878K016	878K018	879K018
Diameter (mm)	1.4	1.6	1.8	1.8
Head Length (mm)	8.0	8.0	8.0	9.5
Super Fine				
Fine				
Medium				Ⓜ
Coarse	C	C	C	
Super Coarse				

TAPER ROUNDED SHOULDER					
Code	845KR018	845KR025	846KR016	847KR016	847KR018
Diameter (mm)	1.8	2.5	1.6	1.6	1.8
Head Length (mm)	4.0	4.0	6.0	8.0	8.0
Super Fine					
Fine		F			
Medium		Ⓜ	Ⓜ	Ⓜ	Ⓜ
Coarse	C	C	C	C	C
Super Coarse					

ORDER INFO:
Please add prefix
70-D
EX: 70-D846012C

ROUND END TAPER												
Code	850012	850014	850016	850018	855012	855014	855018	855025	856010	856012	856014	856016
Diameter (mm)	1.2	1.4	1.6	1.8	1.2	1.4	1.8	2.5	1.0	1.2	1.4	1.6
Head Length (mm)	10.0	10.0	10.0	10.0	6.3	6.3	6.4	7.0	8.0	8.0	8.0	8.0
Super Fine												
Fine			F	F								F
Medium				Ⓜ	Ⓜ	Ⓜ						Ⓜ
Coarse	C	C	C	C		C	C	C	C	C	C	C
Super Coarse												SC

ROUND END TAPER (CONTINUED)					
Code	856018	856025	856L016	856L018	856L020
Diameter (mm)	1.8	2.5	1.6	1.8	2.0
Head Length (mm)	9.0	8.0	9.0	9.0	9.0
Super Fine					
Fine	F				
Medium	Ⓜ				
Coarse	C	C	C	C	C
Super Coarse	SC				

SEPARATION		
Code	889008	889010
Diameter (mm)	0.8	1.0
Head Length (mm)	3.0	4.0
Super Fine		
Fine	F	F
Medium		
Coarse		
Super Coarse		

NEEDLE									
Code	858010	858012	858014	858016	859010	859012	859014	859016	859016
Diameter (mm)	1.0	1.2	1.4	1.6	1.0	1.2	1.4	1.6	1.6
Head Length (mm)	8.0	8.0	8.0	8.0	10.0	10.0	10.0	10.0	10.5
Super Fine			SF				SF	SF	
Fine	F				F		F	F	F
Medium		Ⓜ	Ⓜ	Ⓜ	Ⓜ		Ⓜ	Ⓜ	Ⓜ
Coarse			C			C	C	C	
Super Coarse									

FLAME										
Code	390018	860010	861014	862010	862012	862014	862018	863010	863012	863016
Diameter (mm)	1.8	1.0	1.4	1.0	1.2	1.4	1.8	1.0	1.2	1.6
Head Length (mm)	3.0	4.0	6.0	8.0	8.0	8.0	8.0	10	10	9.5
Super Fine			SF		SF	SF				
Fine	F			F	F	F		F	F	
Medium		Ⓜ			Ⓜ				Ⓜ	Ⓜ
Coarse				C	C	C	C		C	C
Super Coarse										



SF: Super Fine
F: Fine
Ⓜ: Medium
C: Coarse
SC: Super Coarse