

Impression Material Mixing Tip Chart

■ Red requires tip item # 48-778215M

■ Blue requires tip #48-778213M

■ Green requires tip #48-778951M

Dental City Item Number	Dentsply Sirona Part Number	Product Description
49-626600	626600	Reposil Cartridge Light Body Orange 4/Pack
49-678415	678415	Aquasil Smart Wetting LV Teal Fast Set 4 50ml Cartridge Refill 678415
49-678668	678668	Aquasil Ultra+ Smart Wetting XLV Regular Set 4 (50ml) Cartridge Refill
49-678670	678670	Aquasil Ultra+ Smart Wetting XLV Fast Set 4 (50ml) Cartridge Refill
49-678671	678671	Aquasil Ultra+ Smart Wetting LV Regular Set 4 (50ml) Cartridge Refill
49-678672	678672	Aquasil Ultra+ Smart Wetting LV Fast Set 4 (50ml) Cartridge Refill
49-678701	678701	Aquasil Ultra+ XLV Extended Set 4/50ml Cartridges & Tips 678701
49-678771	678771	Aquasil Ultra Refill RS LV 4/Box 678771
49-678772	678772	Aquasil Ultra Refill FS LV 4/Pack 678772
49-678768	678768	Aquasil Ultra Refill RS XLV 4/Box 678768-8
49-678770	678770	Aquasil Ultra Refill FS XLV 4/Box 678770
48-528	619425	Regisil Rigid Refill Cartridge 50ml 619425
48-619510	619510	Regisil PB Cart 4/Pack 619510
48-619500	619500	Regisil 2X Cartridge 4/Pack 619500
44-61E800	61E800	Algin-X Ultra Refill 50ml (4) Fast Set 50ml Cartridge 61E800
44-61E801	61E801	Algin-X Ultra Bulk Fast Set 24/pk 61E801
49-626210	626210	Reposil Cartridge Regular Body Blue 4/Pack
49-626240	626240	Reposil VPS Regular Body Blue 12 Cartridge Bonus Pack 626240
49-626220	626220	Reposil Cartridge Heavy Body Green 4/Pack
49-626250	626250	Reposil VPS Heavy Body 12 Cartridge Bonus Pack
49-626230	626230	Reposil Quixx Putty Cart 4/Pack 626230
49-678056	678056	Aquasil Ultra Cordless Tray Material 4x50ml
49-678250	678250	Aquasil Smart Wetting Monophase RS Purple (4) 50ml 678250
49-678350	678350	Aquasil Smart Wetting Rigid Regular Set 50 ml 4/Box 678350
49-678355	678355	Aquasil Rigid Fast Set 50 ml 4/Box 678355
49-678602	678602	Aquasil Ultra+ Smart Wetting Medium Regular Set 24 (50ml) Cartridge Bulk
49-678673	678673	Aquasil Ultra+ Smart Wetting Medium Regular Set 4 (50ml) Cartridge Refill
49-678601	678601	Aquasil Ultra+ Smart Wetting Medium Fast Set 24 (50ml) Cartridge Bulk
49-678674	678674	Aquasil Ultra+ Smart Wetting Medium Fast Set 4 (50ml) Cartridge Refill
49-678625	678625	Aquasil Ultra+ Smart Wetting Heavy Regular Set 24 (50ml) Cartridge Bulk
49-678675	678675	Aquasil Ultra+ Smart Wetting Heavy Regular Set 4 (50ml) Cartridge Refill
49-678626	678626	Aquasil Ultra+ Smart Wetting Heavy Fast Set 24 (50ml) Cartridge Bulk
49-678676	678676	Aquasil Ultra+ Smart Wetting Heavy Fast Set 4 (50ml) Cartridge Refill
49-678695	678695	Aquasil Ultra+ Smart Wetting Impression Material Intro Kit LV/XLV/Heavy Fast Set
49-678627	678627	Aquasil Ultra+ Smart Wetting Rigid Fast Set 24 (50ml) Cartridge Bulk
49-678678	678678	Aquasil Ultra+ Smart Wetting Rigid Fast Set 4 (50ml) Cartridge Refill
49-678703	678703	Aquasil Ultra+ Heavy Extended Set 4/50ml Cartridges & Tips 678703
49-678773	678773	Aquasil Ultra Refill RS Monophase 4/Box 678773
49-678774	678774	Aquasil Ultra Refill FS Monophase 4/Box 678774
49-678777	678777	Aquasil Ultra Refill RS Rigid 4/Box 678777
49-678778	678778	Aquasil Ultra Refill FS Rigid 4/Box 678778
49-678275	678275	Aquasil Ultra Smart Wetting Tray Heavy-Regular Set Bulk Package 24 x 50mL
49-678775	678775	Aquasil Ultra Refill RS Heavy 4/Box 678775
49-678276	678276	Aquasil Ultra Smart Wetting Tray Heavy- Fast Set Bulk Package 24 x 50mL
49-678776	678776	Aquasil Ultra Refill FS Heavy 4/Box 678776
44-61E810	61E810	Algin-X Ultra Deca Standard Fast Set 2/Pack 61E810
49-678059	678059	Aquasil Ultra Cordless Tray Material 380ml DECA Economy Pack 678059-3
49-678697	678697	Aquasil Ultra+ Smart Wetting Medium Fast Set 2 (380ml) DECA Cartridges Standard
49-678682	678682	Aquasil Ultra+ Smart Wetting Medium Regular Set 2 (380ml) DECA Cartridges Standard
49-678684	678684	Aquasil Ultra+ Smart Wetting Heavy Fast Set 2 (380ml) DECA Cartridges Standard
49-678686	678686	Aquasil Ultra+ Smart Wetting Heavy Regular Set 2 (380ml) DECA Cartridges Standard
49-678688	678688	Aquasil Ultra+ Smart Wetting Rigid Fast Set 2 (380ml) DECA Cartridges Standard
49-678704	678704	Aquasil Ultra+ DECA Heavy Extended Set 2/380ml Cartridges & 20/Tips 678704
49-678788	678788	Aquasil Ultra Deca Standard Fast Set Rigid 2/Pack 678788
49-678782	678782	Aquasil Ultra Deca Monophase 678782
49-678786	678786	Aquasil Ultra Deca Standard Kit Heavy Regular Set 678786
49-678784	678784	Aquasil Ultra Deca Heavy Body Fast Set 678784
49-678797	678797	Aquasil Ultra Deca Monophase Fast Set 678797