




SafeGauze[®] HemoStat[™]

Topical Hemostatic Dressing




Quickly transforms into a clear gel to control bleeding—No suturing required.


HemoStat is a 100% natural, dissolvable hemostatic dressing designed to control bleeding after extractions, periodontal surgery and other procedures that require bleeding control.



Absorbs fluid from oral mucosa and quickly transforms into a clear viscous gel for easy placement.



Creates a gelatinous “plug” that adheres to the wound, creating a tight seal that closes off blood vessels to control bleeding and start the healing process.



Gelatinous plug gradually dissolves into saline and glucose.

See for yourself how well SafeGauze HemoStat works. Order your samples today at Medicom.com or contact your sales representative.



FEATURES	BENEFITS
Transforms into a gel within seconds	Fast and easy placement
Absorbs 400% its weight in fluid	Fills wound and forms a strong clot for effective bleeding control
Adheres to the wound	Stays securely in place—no suturing required
Made from 100% hypo-allergenic cellulose with no chemicals or animal-based products	Safe, natural treatment; helps prevent adverse effects
Natural antibacterial properties	Promotes healing and reduces incidence of post-operative infection
Gelatinous plug gradually dissolves into saline and glucose	Minimizes the risk of infection
Convenient packaging	Saves valuable chair time

What dentists are saying about SafeGauze HemoStat:

“When placed in the wound void, the topical dressing forms a clear gel that expands rapidly as it absorbs blood, creating a plug that stops the bleeding. The dressing then stays in place and does not need to be sutured over.”

“With this product, there’s really no learning curve. You can begin using it right away.”

“Effective, economical and a time saver.”

“Another benefit of SafeGauze HemoStat is that it doesn’t stick to instruments while packing in the socket, which makes it a very quick and easy process. But most importantly, the product controls the bleeding, maintains the clot, and promotes wound healing before eventually dissolving to saline and glucose.”



Canadian code

4900-STATC

USA code

4900-STAT

0.75" x 0.75" - 4 Ply, 20 units per box



Single Use



Not made with natural rubber latex.



Made in China

Rev.01

Like you, we never compromise. Our commitment is to help to make the world safer and healthier by offering consistent, reliable protection. We use carefully selected materials and engineer our products to offer protection you can count on.

As the maker of Safe brands, we take pride in the protection we offer, so you can take pride in yours.

Medicom® Pride in Protection.

www.medicom.com

AMD Medicom Inc.
T. (514) 636-6262 F. (514) 636-6266
cserviceemail@medicom.ca

Toll Free (North America) T. 1 (800) 361-2862 F. 1 (800) 616-4266