

July 13, 2022

Subject: Lead content in Crosstex products

Product Type/Part Number(s): Steam Indicator Strip- 4"/IOcm (SIS-250)

To whom it may concern,

Please be advised that the Crosstex products listed above have been tested by an independent third party. Analytical test results indicate that the lead content in the products is less than 5mg/L. This complies with limits set by US state and federal regulations.

PRODUCT	PART NUMBER(s)	LEAD CONTENT (mg/L)	REGULATORY LEVEL (mg/L)
Steam Indicator Strip- 4"	SIS-250	0.12	5.0

As per 40 CFR 261.24, Steam Indicator Strip- 4" is not considered hazardous waste and can be disposed of in regular trash.

Sincerely,



Regine Celestin
Regulatory Affairs Specialist