

3M Science.
Applied to Life.™

3M™ RelyX™ Unicem and Unicem 2
Self-Adhesive Resin Cements

All-around reliability.




The world's most clinically proven self-adhesive resin cement technology

For simplicity, convenience and reliability, 3M™ RelyX™ Unicem Self-Adhesive Resin Cement is your “go-to” product for a wide range of indications. Thanks to the dual-cure, self-adhesive chemistry, Unicem cements deliver high adhesion, resistance to discoloration, long-term stability and more.

Ideal for:
General cementation, oxide ceramic
restorations, posts and bridges

Results of *The Dental Advisor*
15 year clinical performance study
+++++



Source: *The Dental Advisor* 15-year Clinical Performance Study, March 2017
Rating system: Excellent +++++

Formulated to meet the needs of your practice – and your patients!

- Saves time by eliminating etching, priming and bonding steps
- Virtually no post-op sensitivity
- Strong
- Moisture tolerant
- Endo tips make cementing posts easy and reliable



RelyX™ Unicem 2 Automix Cement Refill – A2 (56846)



RelyX™ Unicem 2 Clicker™ Dispenser Refill – A2 (56875)

Recommended Applications (++ Highly recommended / + Recommended)

Indication	Metal/ Metal Based	Glass Ceramics (including Li. Disilicate) (e.g. e.max®, VITA Mark II, IPS Empress 2)	Oxide Ceramics (e.g. Lava™ Plus, Brux Zir®, Procera®)	Resin Composites (e.g. VITA ENAMIC®)
Inlays/Onlays	++	+	++	+
Crowns/Bridges	++	++	++	+++*
Endodontic Posts	++		++	++
Maryland Bridges	+		+	
On Implant Abutments	++	+	++	++

*crowns only

Ordering Information

Item # Product Information

RelyX™ Unicem 2 Automix Value Packs

Item No.	Product	Contents
56858	Shade Translucent	(3) 8.5g syringes
56849	Shade A2 Universal	(30) mixing tips (regular)
56857	Shade A3 Opaque	(20) mixing tips (wide)
		(15) endo tips; (5) intra-oral tips

RelyX™ Unicem 2 Automix Refills

Item No.	Product	Contents
56848	Shade Translucent	(1) 8.5g syringes
56846	Shade A2 Universal	(10) mixing tips (regular)
		(5) mixing tips (wide)
56847	Shade A3 Opaque	(5) endo tips

RelyX™ Unicem 2 Automix Mix Tips

Item No.	Product
56851	(30) mixing tips (regular)
56918	(15) mixing tips (wide); (15) intra-oral tips
56919	(15) mixing tips (wide); (15) endo tips

RelyX™ Unicem 2 Clicker™ Dispenser Refills

Item No.	Product	Contents
56874	Shade Translucent	(1) 11g Clicker™ Dispenser
56875	Shade A2 Universal	(1) mix pad (70144)
		technique card
56876	Shade A3 Opaque	instructions for use



3M Oral Care
2510 Conway Avenue
St. Paul, MN 55144-1000 USA

Phone 1-800-634-2249
Web 3m.com/dentalcements

3M Canada
Post Office Box 5757
London, Ontario N6A 4T1
Canada

Phone 1-888-363-3685

3M, “3M Science. Applied to Life.”, Clicker and RelyX are trademarks of 3M Company or 3M Deutschland GmbH. Used under license in Canada. All other trademarks are owned by their respective owners.

Please recycle. © 3M 2018.
All rights reserved. 70-2013-0388-3 HB14297 REV 5/2018