

NEW!

# ONE FOR THE BASICS

Venus®

Product Information



[kulzerUS.com/  
try-one](http://kulzerUS.com/try-one)

BISGMA  
& BPA  
FREE!

Giving a hand to oral health.



**KULZER**  
MITSUI CHEMICALS GROUP

## ALWAYS GRAB THE RIGHT SHADE FOR YOUR EVERYDAY CASES.



Venus Diamond/Venus Pearl ONE Shade blends seamlessly into the surrounding dentition, independent of the tooth shade, enabling restorations to simply disappear. The universal composites for everyday restorations offer a new one-shade solution within the Venus family for the majority of everyday cases.



**ONE SHADE FOR ALL –  
THANKS TO THE AMAZING  
ADAPTIVE LIGHT MATCHING.**

Venus Diamond/Venus Pearl ONE Shade determines the shade of the restoration by absorbing the lightwaves that are reflected by the surrounding tooth shade.



**"The One Solution"**  
Start Video:  
[kulzer.com/  
one-effect](https://kulzer.com/one-effect)

## **EVERYDAY RESTORATIONS HAVE NEVER BEEN MORE CONVENIENT.**

The Venus Diamond and Venus Pearl highlights for easy handling:


- Non-slump qualities for easy carving and sculpting
- Stable, minimally adhesive, with an extended working time
- Choose your preferred consistency:
  - The firm, low-stick consistency of Venus Diamond
  - The creamy, low-stick consistency of Venus Pearl



Excavated cavity with pulp capping. Case by: Prof. C.-P. Ernst, Mediplus Mainz, Germany



The restoration matches perfectly into the tooth structure and to both neighboring teeth.




Venus Diamond/Venus Pearl ONE Shade also has a high radiopacity ensuring a reliable diagnosis.  
Source: Kulzer GmbH, SCA, Germany

## LONG-LASTING RESTORATIONS THANKS TO THE UNIQUE TCD-MATRIX.

ONE Shade benefits from the unique combination of low shrinkage stress and high flexural strength of Venus Diamond and Venus Pearl which has been well-documented and tested in clinical studies.

### Venus Diamond & Venus Pearl: Unique combination of low shrinkage stress and high flexural strength.



Source: Kulzer R&D, Wehrheim

Admira Fusion® is a trademark of VOCO GmbH · Tetric EvoCeram® and IPS Empress Direct® are trademarks of Ivoclar Vivadent AG · Estelite® is a trademark of Tokuyama Dental Corporation Inc. · Herculite® is a trademark of Kerr Corporation · Filtek® is a trademark of 3M Company · Spectrum®, TPH® and x-flow® are trademarks of Dentsply DeTrey GmbH · Esthet-X® is a trademark of Dentsply Sirona Inc.

# WANT TO TRY ONE?

We have already packed your ONE Shade kit for you. Order now and try it out for yourself.

kulzerUS.com/  
try-one

Venus®



Art Code	Package	Content
<b>Venus Pearl ONE Shade</b>		
66081836	Venus Pearl Syringe Refill ONE	1 x 3 g
66081837	Venus Pearl PLT Refill ONE	20 x 0.2 g
66081840	Venus Pearl ONE Kit – Syringe	2 x 3 g Syringes of Venus Pearl in shade ONE / 1 x 4 ml bottle of iBOND Universal
66081841	Venus Pearl ONE Kit – PLT	30 x 0.2 g PLT of Venus Pearl in shade ONE / 1 x 4 ml bottle of iBOND Universal
<b>Venus Diamond ONE Shade</b>		
66081838	Venus Diamond Syringe Refill ONE	1 x 4 g
66081839	Venus Diamond PLT Refill ONE	20 x 0.25 g
66081842	Venus Diamond ONE Kit – Syringe	2 x 4 g Syringes of Venus Diamond in shade ONE / 1 x 4 ml bottle of iBOND Universal
66081843	Venus Diamond ONE Kit – PLT	30 x 0.25 g PLT of Venus Diamond in shade ONE / 1 x 4 ml bottle of iBOND Universal

## Distributed by:

Kulzer, LLC  
4315 S. Lafayette Blvd.,  
South Bend, IN 46614  
(800) 431-1785

## Manufactured by:

Kulzer GmbH  
Leipziger Straße 2  
63450 Hanau  
Germany  
info.dent@kulzer-dental.com

kulzerUS.com

© 2020 Kulzer GmbH. All Rights Reserved.

ONE\_BRO\_108/2020