



VersaSure™ Cleaner Disinfectant Wipes

Kills viruses* and bacteria in 2 minutes



Versatility to use everywhere with the assurance of broad spectrum disinfection

Efficacy

- ▶ Kills bacteria and fungus, including MRSA and TB, in 2 min
- ▶ Kills viruses, including mumps and measles, in 30 sec†

Aesthetics

- ▶ Alcohol-free
- ▶ Low odor
- ▶ Low residue

Compatibility

- ▶ Good compatibility on surfaces commonly found in the healthcare setting



	VersaSure™	PDI® Sani-Cloth® AF3	PDI® Super Sani-Cloth®	CaviWipes
Bacteria and virus* kill times‡	30 sec—2 min	3 min	2 min	3 min
Total kill claims‡	44	44	29	14
Toxicity Signal Word	Caution	Caution	Warning	Caution

* See virus claims table on product †Rotavirus is 2 minutes

‡ EPA-registered disinfection claims listed on product labels as of April 2019. For a full list of all EPA-registered kill claims and contact times, refer to product Master Label.

Clorox Healthcare® VersaSure™ Kill Claims: 44 pathogens in 2 minutes or less

Bacteria - 2 minutes
<i>Bordetella pertussis</i>
<i>Burkholderia cepacia</i>
<i>Campylobacter jejuni</i>
Carbapenem-resistant <i>Escherichia coli</i>
Carbapenem-resistant <i>Klebsiella pneumoniae</i>
Community-associated methicillin-resistant <i>Staphylococcus aureus</i> (CA-MRSA)*
<i>Enterobacter aerogenes</i>
<i>Escherichia coli</i> (E.coli)
<i>Escherichia coli</i> O157:H7
<i>Klebsiella pneumoniae</i>
<i>Legionella pneumophila</i>
Linezolid-resistant <i>Staphylococcus aureus</i> (LRSA)
Methicillin-resistant <i>Staphylococcus aureus</i> (MRSA)
Multidrug-resistant [MDR] <i>Acinetobacter baumannii</i>
Multidrug-resistant [MDR] <i>Enterococcus faecium</i>
Multidrug-resistant <i>Klebsiella pneumoniae</i>
<i>Mycobacterium bovis-BCG</i> (TB)
Penicillin-resistant <i>Streptococcus pneumoniae</i>
<i>Pseudomonas aeruginosa</i> (Pseudomonas)
<i>Salmonella enterica</i> (Salmonella)
<i>Serratia marcescens</i>
<i>Shigella dysenteriae</i>
<i>Staphylococcus aureus</i>
<i>Staphylococcus epidermidis</i> (Coagulase-negative staphylococci)
<i>Streptococcus pyogenes</i>
Vancomycin-intermediate-resistant <i>Staphylococcus aureus</i> (VISA)
Vancomycin-resistant <i>Enterococcus faecalis</i> (VRE)
Vancomycin-resistant <i>Staphylococcus aureus</i> (VRSA)
Viruses - 30 seconds
Avian influenza Type A Virus (H5N1)
Bovine Viral Diarrhea Virus
Cytomegalovirus
Duck Hepatitis B virus as surrogate for Human Hepatitis B
Herpes Simplex Virus Type 1
Herpes Simplex Virus Type 2
Human Immunodeficiency Virus Type 1
Influenza A (H1N1) Virus
Influenza A Virus
Measles Virus
MERS (Middle East Respiratory Syndrome)-Associated Coronavirus (MERS-CoV)
Mumps virus
Respiratory Syncytial Virus (RSV)
Viruses - 2 minute
Rotavirus
Fungi - 2 minute
<i>Candida albicans</i>
TOTAL: 44 CLAIMS

*US 300 and US 400 strains

Features patented technology that enhances quat activity on surface to provide faster kill times without alcohol

Designed for Use on Common Healthcare Surfaces

- ▶ Common plastics
- ▶ Stainless steel
- ▶ Glass
- ▶ High-end finishes (Corian™, marble)

Ideal for Use in Patient Areas

- ▶ Little residue
- ▶ Excellent cleaning power
- ▶ Stays wet for full 2 minutes
- ▶ Minimal odor

VersaSure™ Cleaner Disinfectant Wipes FOR HARD SURFACES

Multipurpose Wipes

Designed for: *Equipment and hard-surface disinfecting*

Item No. 31757, 85 ct. canister, 6/case (6.75" x 8")

Item No. 31758, 150 ct. canister, 6/case (6" x 5")

Item No. 31760, 30 ct. pouch, 24/case (6.75" x 8")



Terminal Wipes

Designed for: *Everyday cleaning and disinfecting, daily patient room cleaning*

Item No. 31759, 110 ct. bucket, 2/case (12" x 12")

Item No. 31761, 110 ct. refill pouch, 2/case (12" x 12")



HEALTHCARE™

FROM CLOROX PRO™

For product resources and implementation tools, contact your Clorox sales representative or Call: 1-800-234-7700 Visit: www.CloroxHealthcare.com Email: healthcare@clorox.com

© 2020 Clorox Professional Products Company, 1221 Broadway, Oakland, CA 94612.